


Introduction: Despite notorious evolution, nebulized therapies continue to face important unmet needs, particularly when considering the following three important concepts:



1. **Uncontrolled particle deposition** remains a limitation for effective targeted drug delivery into the lungs. The **breathing manoeuvre** has a significant impact on deposition, especially when smaller airways are targeted. Therefore, drug delivery is frequently **uncontrolled** or not tailored to the targeted deposition areas due to a large variability in the inhalation manoeuvre of nebulizer users. This leads to **lack of efficacy** in pulmonary drug deposition and consequently to **higher risks and costs**.



2. High risks in clinical development are associated to the inability to monitor **high influential parameters** on drug deposition, adding difficulty and complexity to clinical trials. Inefficient management of clinical/device data also poses a challenge to effective study monitoring. In addition, the considerable results variability often require **costly study repetitions** for dose finding and further safety evaluations, thereby increasing overall expenses.



3. The patient's needs include:

- improving the nebulizer's **usability** by optimizing its assembly, use and cleaning process;
- a device that provides **feedback and guidance** for the patient to achieve an **optimal inhalation manoeuvre**;
- an **alleviation in the disease burden**, by simplifying the interaction with the nebulizer and facilitating therapy adherence and compliance,

Overall, it's crucial to find a **balance** between treatment effectiveness and patient engagement.

Results: Focused on these challenges, Pulmotree developed Kolibri, a connected mesh nebulizer, whose technology is hereby described.



User-specific Customization
for medication, patient group and market strategy



Haptic and digital Inhalation feedback
for optimized lung deposition



Intelligent Administration System (iNAS)
for breath-triggered aerosol generation and targeted drug deposition




Automatic data transfer to users, physicians and CROs, with the Pulmotree Connect



Key- & Lock-System
for product protection and easy handling

Kolibri as a nebulizer platform can be **customized** according to each drug-device combination, due to a **tailored aerosol generator** that produces optimal aerosol particles.

Additionally, an **unique inhalation technology** guides the patient through **haptic and digital feedback** to achieve the ideal inspiration manoeuvre. For optimized lung deposition, the breath-triggered aerosol generation and the **iNAS technology** only administer aerosol in the therapeutically useful phase. By analysing the patient's behaviour, iNAS trains the patient to perform the desired inhalation manoeuvre.



Remote monitoring is possible due to **automatic data transfer** on inhalation, adherence and device status, while using a centralized interface platform and a mobile app.

A **key- & lock-system** protects the product on the market and also protects the patient of interchangeability with other medications.

For an enhanced **patient's experience and usability**, Kolibri is a hand-held small portable device of simple assembly and cleaning, with a modern design that alleviates the treatment's onus.

Conclusion: Kolibri's technology presents solutions to unanswered challenges and aims to drive progress in nebulized therapies.

Furthermore, by providing **new patient-side parameters** generated during nebulization, important **inputs** would be available for **lung deposition models**. These future developments can mark a significant advancement in personalized respiratory medicine, fostering further knowledge and pushing evolution in inhaled therapies.

